

## **ARTONE Communication Solutions**

## Artone MIC

ARTONE MIC Utilize Bluetooth headphones or earbuds as an personalized hearing amplifier. Simply pair Bluetooth headphones with Artone MIC and enjoy amplified sound like never before. Adjust the sensitivity based on surroundings, whether it's meetings, restaurants, or car rides. In any situation where hearing becomes challenging, just hand over the Artone MIC to the speaker and get ready to catch every word with clarity.

- Crystal-Clear Phone Calls: Enjoy hands-free conversations with exceptional sound quality.
  Artone-3 MAX eliminates background noise, ensuring every word comes through loud and clear.
  (Keywords: phone calls, hands-free, hearing aid, cochlear implant)
- Stream What You Love: Listen to music, podcasts, and audiobooks directly from your phone or tablet. Stream audio, music streaming, audiobooks, podcasts, without any interference.
- Turn the TV into a Personal Oasis: Ditch the bulky headphones and experience rich, clear TV sound streamed directly to a hearing aid or Cochlear implant, minimize the impact of hearing loss on the quality of life.
- Universal Compatibility: Artone-3 MAX works seamlessly with virtually any Bluetooth device and most T-coil equipped hearing aids and implants
- No feedback, no squeaking, no metallic voice, no oclusion "barrel" effect.
- Cristal clear sound
- Bluetooth V4.1 and supports profiles A2DP/ HS/ HF/ AVRCP
- Compatibility: Smartphone | Bluetooth Hearing Aids | Cochlear implants | Earphones/ Buds
- electively captures the sound you need, transmitting it directly to your Bluetooth headphones when clipped to the speaker's shirt.
- Focus on the speech band 8kHz

• Dimmensions: 53mm x 18.3mm x 17mm | Weight: 9.5gr

Range: 10 - 20 meter in/ ourdoors

Power: Rechargeable battery

Battery last time: Up to 7 hours

Charging time: 90minCharging voltage: 5V DC

Wireless Version: V 4.1

