

ARTONE Communication Solutions

Artone MS - Bluetooth Neckloop

Artone MS is an advanced version of the Wireless Bluetooth Loopset that supports two-link signal connection and A2DP transmission. Thanks to these features, the ARTONE MS bluetooth neckloop enables hearing aid user to answer calls on their mobile phone while watching TV simply by pressing one button.

- ARTONE MS is equipped with a neckloop that allows to connect the device to any T-Coil hearing aid. Thereby, the ARTONE MS turns the existing hearing aid into a wireless communicator that connects to any Bluetooth-enabled audio source such as mobile phone, music player, PC. To use the ARTONE MS with non-Bluetooth devices, an additional transmitter (ARTONE TVB or ARTONE TVB+) must be used.
- Using the ARTONE MS, the hearing aid user can communicate over their mobile phone using their hearing aid as handsfree or listen to music using their hearing aid as wireless headphones.
- Crystal clear amplified sound of voice or music is transmitted from the audio source to the hearing aid through the ARTONE MS without any interference.
- ARTONE MS will alert the hearing aid user about an incoming call while interrupting the TV sound
 - 2 Bluetooth device connectivity.
 - Stereo frequency range - A2DP
 - Totally wireless, hands-free
 - Clear amplified communication.
 - Compatible with any phone with Bluetooth® technology.
 - Compatible with any hearing aid equipped with T-coil option
 - Rechargeable Long lasting Battery
 - Voice dial options.
 - Works MP3 or any other music device that is Bluetooth® compatible.
(compatible with Artone TVB and Artone TVB+ – for SKYPE, TV, Radio, Ipod...)
- Light and elegant style. ARTONE MS will then let the hearing aid user switch to conversation mode with a press of one button. When the call is finished, the signal of the TV will get back to the hearing aid.

