

# PHONICOMM<sup>®</sup> | JAZZ hearing amplifier

Phonicomm Jazz is a Thin Tube self-fitting OTC hearing amplifier providing a do-it-yourself hearing solution which can be fitted and adjusted to personal preferences. Developed with the understanding and knowledge of mild to moderate hearing losses, Jazz hearing amplifiers prioritize speech, customizable settings, ease of use, comfort, convenience, and Phonicomm signature sound quality. Packed in an attractive housing, Jazz offers a tailor-made hearing solution for people facing problems with beginning hearing loss and may be the first step to conventional hearing aids.

## Features:

- Monaural
- Advanced speech enhancement
- Noise Reduction
- Active feedback cancelation
- Easy setup and personalization from the Earfit App
- Connectable to any Bluetooth smartphone, laptops, tablet, smart TV, etc.)
- Easy to Use
- Rechargeable
- Charging time 2 hours



## Tech Spec:

OSPL 90	125dB + 3dB
HFA OSPL 90	111dB ± 4dB
HFA Full-On Gain	400 dB ± 5dB
Frequency Respons	350 Hz - 4500 Hz
Total Harmonic Distortion	<4% + 3%
	500 Hz - 2.0%
	800 Hz - 0.3%
	1600 Hz - 0.6%
Equivalent Input Noise	<23dB + 3dB
Channels	16
Frequency bands	64
Noise Reduction	Yes
Feedback Canselation	Yes
Bloetooth	Yes
App	Yes